


# Scanning Instructions

MOVEMENT HEALTH PLATFORM



## CREATE OR FIND YOUR PROFILE

### First Time?

1. Open the Sparta Scan application
2. Click Add Users icon 
3. Select Team in the drop-down
4. Enter all information available
5. Click Create User and Continue

### Returner?

1. Login to the Sparta Scan App
2. Select Team from dropdown list
3. Find and select profile
4. Select Scan type
5. Click Start Scan

Start Scan

## SCAN PROTOCOLS

### Progressive Balance Scan

1. Find and select a user, click on Balance, and click Start Scan
2. Follow the prompts on the screen:
  - a. Step on the plate and stand still
  - b. When prompted, shift your weight to your RIGHT until the dot on the screen is within the highlighted area, and hold
  - c. At the tone, return your weight to the center and stand still
  - d. When prompted, shift your weight to your LEFT until the dot on the screen is within the highlighted area, and hold
  - e. At the tone, return your weight to the center and stand still
  - f. When prompted, lift your left leg to balance on the RIGHT for 8-10 seconds
  - g. Return your left foot to the plate and stand still
  - h. When prompted, repeat the above steps on the LEFT side
  - i. Continue following prompts for a total of two balance trials on each side
3. If a Balance task cannot be completed successfully, click "Skip Segment" or "Skip All" to complete the scan

### Balance Scan

1. Find and select a user, click on Balance, and click Start Scan
2. Follow the prompts on the screen:
  - a. Step on the plate and stand still
  - b. When prompted, lift your left leg to balance on the RIGHT for 8-10 seconds
  - c. Return your left foot to the plate and stand still
  - d. When prompted, repeat the above steps on the LEFT side
  - e. Continue following prompts for a total of two balance trials on each side