


# Welcome to Your New Scan Kit



Your Scan Kit includes the following:

- iPad
  - Force Plate
- **A** U-Shape Adapter
  - **B** 2-to-1 Adapter
  - **C** Charge Cable (USB-C)
  - **D** Data Cable (USB-C to USB-B)
  - **E** Power Adapter



Follow the instructions below to get up and running:

## FULL SETUP (Recommended)

1. Plug the male end of the U-Shape Adapter **A** into the iPad, then attach the 2-to-1 Adapter **B** into the female port of the U-Shape Adapter so both adapters wrap behind the iPad.
2. Connect the blocky USB-B connector of the Data Cable **D** into the port on the underside of the Force Plate, then connect the other end to the 2-to-1 Adapter behind the iPad.
3. Connect the Charge Cable **C** to the Power Adapter **E** and the 2-to-1 Adapter, and make sure the Power Adapter is plugged into the wall.
4. Turn on the iPad (using the button opposite the cables) and check for the Sparta Scan app on the home screen. If the app isn't installed, download it from the App Store.
5. Open the Settings app, scroll down, and tap "Sparta Scan" in the left panel. Then, turn on the Force Plate Driver on the right.
6. Open Sparta Scan, and sign in to start scanning.

## QUICK SETUP

Plug the iPad directly to the Force Plate using the Data Cable, then follow steps 4-6 above to start scanning. Please note, you'll need to connect the iPad to power to charge.



**Next:** create an account at [academy.spartascience.com](https://academy.spartascience.com) to get started using the Sparta Movement Health Platform.